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
3) Espectro Emisión HLED a 625 nm (-+10%)

Longitud de onda

Low-level light therapy for androgenetic alopecia: a 24-week, randomized, double-blind, sham device-controlled multicenter trial.

[Kim H¹](#), [Choi JW](#), [Kim JY](#), [Shin JW](#), [Lee SJ](#), [Huh CH](#).

+ Author information

Abstract

BACKGROUND: Androgenetic alopecia (AGA) is a common disorder affecting men and women. Finasteride and minoxidil are well-known, effective treatment methods, but patients who exhibit a poor response to these methods have no additional adequate treatment modalities.

OBJECTIVE: To evaluate the efficacy and safety of a low-level light therapy (LLLT) device for the treatment of AGA.

METHODS: This study was designed as a 24-week, randomized, double-blind, sham device-controlled trial. Forty subjects with AGA were enrolled and scheduled to receive treatment with a helmet-type, home-use LLLT device emitting wavelengths of 630, 650, and 660 nm or a sham device for 18 minutes daily. Investigator and subject performed phototrichogram assessment (hair density and thickness) and global assessment of hair regrowth for evaluation.


RESULTS: After 24 weeks of treatment, the LLLT group showed significantly greater hair density than the sham device group. Mean hair diameter improved statistically significantly more in the LLLT group than in the sham device group. Investigator global assessment showed a significant difference between the two groups, but that of the subject did not. No serious adverse reactions were detected.

CONCLUSION: LLLT could be an effective treatment for AGA.

© 2013 by the American Society for Dermatologic Surgery, Inc. Published by Wiley Periodicals, Inc.

3) Espectro Emisión HLED a 625 nm (-+10%)

Longitud de onda

3) Espectro Emisión HLED a 625 nm (-+10%)

Longitud de onda

LED (660 nm) and laser (670 nm) use on skin flap viability: angiogenesis and mast cells on transition line.

Nishioka MA¹, Pinfildi CE, Sheliga TR, Arias VE, Gomes HC, Ferreira LM.

⊕ Author information

Abstract


Skin flap procedures are commonly used in plastic surgery. Failures can follow, leading to the necrosis of the flap. Therefore, many studies use LLLT to improve flap viability. Currently, the LED has been introduced as an alternative to LLLT. The objective of this study was to evaluate the effect of LLLT and LED on the viability of random skin flaps in rats. Forty-eight rats were divided into four groups, and a random skin flap (10 × 4 cm) was performed in all animals. Group 1 was the sham group; group 2 was submitted to LLLT 660 nm, 0.14 J; group 3 with LED 630 nm, 2.49 J, and group 4 with LLLT 660 nm, with 2.49 J. Irradiation was applied after surgery and repeated on the four subsequent days. On the 7th postoperative day, the percentage of flap necrosis was calculated and skin samples were collected from the viable area and from the transition line of the flap to evaluate blood vessels and mast cells. The percentage of necrosis was significantly lower in groups 3 and 4 compared to groups 1 and 2. Concerning blood vessels and mast cell numbers, only the animals in group 3 showed significant increase compared to group 1 in the skin sample of the transition line. LED and LLLT with the same total energies were effective in increasing viability of random skin flaps. LED was more effective in increasing the number of mast cells and blood vessels in the transition line of random skin flaps.

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